


Product Information	
Catalog No.	PA13428
Product name	Anti-LAMC1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Laminin-7 subunit gamma, Laminin-1 subunit gamma, Laminin subunit gamma-1, Laminin-10 subunit gamma, Laminin-2 subunit gamma, S-laminin subunit gamma, LAMC1, Laminin-11 subunit gamma, Laminin-3 subunit gamma, S-LAM gamma, Laminin-6 subunit gamma, Laminin-9 subunit gamma, Laminin-4 subunit gamma, Laminin-8 subunit gamma, LAMB2, Laminin B2 chain
Immunogen	E. coli - derived recombinant Human LAMC1 (Gly70-Lys287).
Accession	P11047
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-LAMC1 polyclonal antibody at 1ug/mL on recombinant LAMC1 (Catalog No. HY379032).

Western blot