

Product Information	
Catalog No.	PA13452
Product name	Anti-EFNA1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	EPLG1, EFNA1, LERK1, Immediate early response protein B61, EPH-related receptor tyrosine kinase ligand 1, TNFAIP4, TNF alpha-induced protein 4, LERK-1, Ephrin-A1, Tumor necrosis factor alpha-induced protein 4
Immunogen	E. coli - derived recombinant Human EFNA1 (Asp19-Ser182).
Accession	P20827
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.