


Product Information	
Catalog No.	PA13458
Product name	Anti-CD49f/ITGA6 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CD49 antigen-like family member F, Integrin alpha-6, ITGA6, VLA-6, CD49f, Alpha6p
Immunogen	E. coli - derived recombinant Human CD49f/ITGA6 (Val463-Asp1011).
Accession	P23229
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-ITGA6 polyclonal antibody at 1ug/mL on various samples.

- Lane 1: Mouse spleen lysate
- Lane 2: Mouse lung lysate
- Lane 3: Mouse skin lysate
- Lane 4: Rat kidney lysate
- Lane 5: Rat lung lysate
- Lane 6: Rat skin lysate