

Product Information	
Catalog No.	PA13480
Product name	Anti-Human TGFBR2/TGFR-2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	TGFR-2, TbetaR-II, TGF-beta type II receptor, TGFBR2, TGF-beta receptor type-2, TGF-beta receptor type II, Transforming growth factor-beta receptor type II
Immunogen	E. coli - derived recombinant Human TGFBR2/TGFR-2 (Gln49-Glu151).
Accession	P37173
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.