

Product Information	
Catalog No.	PA13485
Product name	Anti-Human CDKN2A/p16INK4a Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	p16-INK4, Multiple tumor suppressor 1, MTS-1, MTS1, Cyclin-dependent kinase 4 inhibitor A, p16-INK4a, Cyclin-dependent kinase inhibitor 2A, CDKN2A, p16INK4A, CDKN2, CDK4I
Immunogen	E. coli - derived recombinant Human CDKN2A/p16INK4a (Met1-Asp156).
Accession	P42771
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.