

Product Information	
Catalog No.	PA13541
Product name	Anti-Human TCIRG1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	TCIRG1, Vacuolar proton translocating ATPase 116 kDa subunit a isoform 3, V-ATPase 116 kDa subunit a3, V-type proton ATPase 116 kDa subunit a3, ATP6V0A3, ATP6N1C, Osteoclastic proton pump 116 kDa subunit, TIRC7, T-cell immune response cDNA7 protein, T-cell immune regulator 1, OC-116 kDa, OC116
Immunogen	E. coli - derived recombinant Human TCIRG1 (Met1-Glu385).
Accession	Q13488
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.