

Product Information	
Catalog No.	PA13554
Product name	Anti-Human RABEP1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	RABPT5, RABEP1, Rabaptin-4, RABPT5A, Rab GTPase-binding effector protein 1, Renal carcinoma antigen NY-REN-17, RAB5EP, Rabaptin-5alpha, Rabaptin-5
Immunogen	E. coli - derived recombinant Human RABEP1 (Ser464-Glu552).
Accession	Q15276
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.