

Product Information	
Catalog No.	PA13618
Product name	Anti-Human AKAP1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	PRKA1, S-AKAP84, AKAP149, AKAP 149, Protein kinase A-anchoring protein 1, D-AKAP-1, Dual specificity A-kinase-anchoring protein 1, AKAP1, Spermatid A-kinase anchor protein 84, A-kinase anchor protein 149 kDa, A-kinase anchor protein 1, mitochondrial
Immunogen	E. coli - derived recombinant Human AKAP1 (Asp606-Leu903).
Accession	Q92667
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.