

Product Information	
Catalog No.	PA13625
Product name	Anti-Human Pinin/PNN/DRSP Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	PNN, MEMA, SR-like protein, Nuclear protein SDK3, Domain-rich serine protein, DRSP, 140 kDa nuclear and cell adhesion-related phosphoprotein, Pinin, Desmosome-associated protein, Melanoma metastasis clone A protein, DRS protein, DRS
Immunogen	E. coli - derived recombinant Human Pinin/PNN/DRSP (Arg123-Arg315).
Accession	Q9H307
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.