


Product Information	
Catalog No.	PA13640
Product name	Anti-Human EXOG Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	ENDOGL2, Endo G-like 1, EXOG, Nuclease EXOG, mitochondrial, ENGL, Endonuclease G-like 1, ENDOGL1
Immunogen	E. coli - derived recombinant Human EXOG (Thr78-Thr298).
Accession	Q9Y2C4
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-EXOG polyclonal antibody at 1µg/ml on various samples.
Lane 1: MCF-7 cell lysate

Western blot