

Product Information	
Catalog No.	PA13650
Product name	Anti-Human SLC32A1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	GABA and glycine transporter, Vesicular inhibitory amino acid transporter, VGAT, Solute carrier family 32 member 1, Vesicular GABA transporter, hVIAAT, SLC32A1, VIAAT
Immunogen	E. coli - derived recombinant Human SLC32A1 (Arg30-Lys115).
Accession	Q9H598
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.