


Product Information	
Catalog No.	PA13742
Product name	Anti-Mouse CD324/CDH1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Cadherin-1 (ARC-1) (Epithelial cadherin) (E-cadherin) (Uvomorulin) (CD antigen CD324) [Cleaved into: E-Cad/CTF1, E-Cad/CTF2, E-Cad/CTF3]
Immunogen	E. coli - derived recombinant Mouse CD324/CDH1 (Ser164-Asn371).
Accession	P09803
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-CDH1 polyclonal antibody at 1µg/ml on various samples.
Lane 1: Mouse testis lysate

Western blot