


Product Information	
Catalog No.	PA13814
Product name	Anti-Mouse TrkB /NTRK2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	BDNF/NT-3 growth factors receptor (EC 2.7.10.1) (GP145-TrkB/GP95-TrkB) (Trk-B) (Neurotrophic tyrosine kinase receptor type 2) (TrkB tyrosine kinase)
Immunogen	E. coli - derived recombinant Mouse TrkB /NTRK2 (Cys32-His377).
Accession	P15209
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-NTRK2 polyclonal antibody at 1µg/ml on various samples.
Lane 1: Mouse brain lysate
Lane 2: Mouse cerebellum lysate
Lane 3: Rat brain lysate
Lane 4: Rat cerebellum lysate