

Product Information	
Catalog No.	PA13855
Product name	Anti-Mouse VSIR/VISTA Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	V-type immunoglobulin domain-containing suppressor of T-cell activation (Platelet receptor Gi24) (V-set domain-containing immunoregulatory receptor) (V-set immunoregulatory receptor)
Immunogen	E. coli - derived recombinant Mouse VSIR/VISTA (Ala32-Asn181).
Accession	Q9D659
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.