


Product Information	
Catalog No.	PA13888
Product name	Anti-Mouse ST6GALNAC2/SIAT7-B Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Alpha-N-acetylgalactosaminide alpha-2, 6-sialyltransferase 2 (EC 2.4.3.3) (Gal-beta-1, 3-GalNAc alpha-2, 6-sialyltransferase) (GalNAc alpha-2, 6-sialyltransferase II) (ST6GalNAc II) (ST6GalNAcII) (Sialyltransferase 7B) (SIAT7-B)
Immunogen	E. coli - derived recombinant Mouse ST6GALNAC2/SIAT7-B (Thr97-Arg373).
Accession	P70277
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-ST6GALNAC2 polyclonal antibody at 1ug/mL on various samples.
Lane 1: Mouse liver lysate

Western blot