

Product Information	
Catalog No.	PA13914
Product name	Anti-HMPV F/Fusion glycoprotein F0 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human metapneumovirus (strain CAN97-83) (HMPV)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Fusion glycoprotein F0, Protein F, Fusion glycoprotein F2, Fusion glycoprotein F1, F
Immunogen	E. coli - derived recombinant HMPV F/Fusion glycoprotein F0 (Phe103-Thr489).
Accession	Q6WB98
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.