

Product Information	
Catalog No.	PA13956
Product name	Anti-Human GRIN1/NMDAR1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	N-methyl-D-aspartate receptor subunit NR1, GluN1, Glutamate [NMDA] receptor subunit zeta-1, Glutamate receptor ionotropic, NMDA 1, NMD-R1, GRIN1, NMDAR1
Immunogen	E. coli - derived recombinant Human GRIN1/NMDAR1 (Asp23-His371).
Accession	Q05586
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.