

Product Information	
Catalog No.	PA13971
Product name	Anti-Escherichia coli pcnB/Poly(A) polymerase I Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Escherichia coli (strain K12)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Poly(A) polymerase I, PAP I, 2.7.7.19, Plasmid copy number protein, pcnB, b0143, JW5808, Escherichia coli (strain K12)
Immunogen	E. coli - derived recombinant Escherichia coli pcnB/Poly(A) polymerase I (Leu13-Ala465).
Accession	P0ABF1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.