

Product Information	
Catalog No.	PA14045
Product name	Anti-Human PDSS2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	All-trans-decaprenyl-diphosphate synthase subunit 2, Solanesyl-diphosphate synthase subunit 2, C6orf210, PDSS2, DLP1, Decaprenyl-diphosphate synthase subunit 2, Candidate tumor suppressor protein, All trans-polyprenyl-diphosphate synthase PDSS2, Decaprenyl pyrophosphate synthase subunit 2
Immunogen	E. coli - derived recombinant Human PDSS2 (Ser143-Ser399).
Accession	Q86YH6
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.