


Product Information	
Catalog No.	PA14083
Product name	Anti-Human TMED9 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Glycoprotein 25L2, GP25L2, p24alpha2, GMP25, p24 family protein alpha-2, p25, Transmembrane emp24 domain-containing protein 9, TMED9
Immunogen	E. coli - derived recombinant Human TMED9 (Leu38-Gly146).
Accession	Q9BVK6
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-TMED9 polyclonal antibody at 1ug/mL on various samples.

Lane 1: HepG2 cell lysate

Lane 2: THP-1 cell lysate

Lane 3: HepG2 cell lysate

Lane 4: MCF-7 cell lysate