


Product Information	
Catalog No.	PA14149
Product name	Anti-Mouse GFAP Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	GFAP, Glial fibrillary acidic protein
Immunogen	E. coli - derived recombinant Mouse GFAP (Thr107-Glu210).
Accession	P03995
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-GFAP polyclonal antibody at 1ug/mL on various samples.

Lane 1: Mouse brain lysate

Lane 2: Mouse cerebellum lysate

Lane 3: Rat brain lysate

Lane 4: Rat cerebellum lysate