


Product Information	
Catalog No.	PA14195
Product name	Anti-Mouse FABP5 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Psoriasis-associated fatty acid-binding protein homolog, FABP5, E-FABP, Epidermal-type fatty acid-binding protein, Fatty acid-binding protein, epidermal, Fatty acid-binding protein 5, PA-FABP
Immunogen	E. coli - derived recombinant Mouse FABP5 (Met1-Gln135).
Accession	Q05816
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-FABP5 polyclonal antibody at 1ug/mL on various samples.
Lane 1: Mouse brain lysate
Lane 2: Mouse skin lysate
Lane 3: Rat brain lysate
Lane 4: Rat liver lysate