


Product Information	
Catalog No.	PA14198
Product name	Anti-Mouse CD244 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	SLAM family member 4, SLAMF4, h2B4, Signaling lymphocytic activation molecule 4, Natural killer cell receptor 2B4, NKR2B4, NAIL, NK cell type I receptor protein 2B4, CD244, NK cell activation-inducing ligand, 2B4
Immunogen	E. coli - derived recombinant Mouse CD244 (Asp24-Thr212).
Accession	Q07763
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-CD244 polyclonal antibody at 1ug/mL on various samples.
Lane 1: MCF-7 cell lysate

Western blot