


Product Information	
Catalog No.	PA14338
Product name	Anti-Human IGFBP1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	IGFBP-1, IGFBP1, Insulin-like growth factor-binding protein 1, Placental protein 12, IBP1, IGF-binding protein 1, IBP-1, PP12
Immunogen	E. coli - derived recombinant Human IGFBP1 (Ala26-Asn259).
Accession	P08833
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-IGFBP1 polyclonal antibody at 1ug/mL on various samples.
Lane 1: HepG2 cell lysate

Western blot