

Product Information	
Catalog No.	PA14425
Product name	Mouse Anti-Romiplostim Polyclonal Antibody
Applications	ELISA, WB
Species reactivity	Human
Host species	Mouse
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Thrombopoietin, C-mpl ligand, ML, Megakaryocyte colony-stimulating factor, Megakaryocyte growth and development factor, MGDF, Myeloproliferative leukemia virus oncogene ligand, THPO, MGDF
Immunogen	Romiplostim
Accession	P40225
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Alternative Names	Romiplate, Romiplostim, Nplate
Product Usage Information	ELISA:1:5000-1:20000;WB:1:500-1:2000
Note	For research use only.