


Product Information	
Catalog No.	PA14428
Product name	Anti-Human MT-ATP8 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	ATP synthase protein 8, A6L, F-ATPase subunit 8, MT-ATP8, ATP8, ATPASE8, MTATP8
Immunogen	E. coli - derived recombinant Human MT-ATP8 (Gln25-Ser68).
Accession	P03928
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-MT-ATP8 polyclonal antibody at 1ug/mL on various samples.

Lane 1: HepG2 cell lysate

Lane 2: Lncap cell lysate

Lane 3: A549 cell lysate

Lane 4: Hela cell lysate