

Product Information	
Catalog No.	PA14430
Product name	Anti-EBV/HHV4 EBNA1/BKRF1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Epstein-Barr nuclear antigen 1, EBNA-1, EBV nuclear antigen 1, EBNA1, BKRF1
Immunogen	E. coli - derived recombinant EBV/HHV4 EBNA1/BKRF1 (His468-Pro607).
Accession	P03211
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.