

Product Information	
Catalog No.	PA14433
Product name	Anti-Imdevimab Polyclonal Antibody
Applications	ELISA
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Rabbit polyclonal to REGN-10987.
Immunogen	Imdevimab
Accession	CAS: 2415933-40-1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	The product can be stored for 2 - 3 weeks at 2 to 8°C or for up to 12 months at -20°C. Avoid repeated freezing and thawing cycles.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Background

Casirivimab (REGN-10933) / Imdevimab (REGN-10987) (Ronapreve™; REGEN-COV™) is a co-packaged combination of two neutralizing immunoglobulin gamma 1 (IgG1) human monoclonal antibodies (casirivimab and imdevimab) against the spike protein of the severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), the causative agent of coronavirus disease 2019 (COVID-19).