

Product Information	
Catalog No.	PA14490
Product name	Anti-Baker's yeast SSB1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Ribosome-associated molecular chaperone SSB1, 3.6.4.10, Cold-inducible protein YG101, Heat shock protein SSB1, Hsp70 chaperone Ssb, SSB1, YG101, YDL229W
Immunogen	E. coli - derived recombinant Baker's yeast SSB1 (Met1-Arg613).
Accession	P11484
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.