

Product Information	
Catalog No.	PA14550
Product name	Anti-Human UTY Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Histone demethylase UTY, 1.14.11.68, Ubiquitously-transcribed TPR protein on the Y chromosome, Ubiquitously-transcribed Y chromosome tetratricopeptide repeat protein, [histone H3]-trimethyl-L-lysine(27) demethylase UTY, UTY, KDM6C
Immunogen	E. coli - derived recombinant Human UTY (Lys1042-Trp1205).
Accession	O14607
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.