


Product Information	
Catalog No.	PA14580
Product name	Anti-Human POLR3G Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	DNA-directed RNA polymerase III subunit RPC7, RNA polymerase III subunit C7, DNA-directed RNA polymerase III subunit G, RNA polymerase III 32 kDa alpha subunit, RPC32-alpha, RNA polymerase III 32 kDa subunit, RPC32, POLR3G
Immunogen	E. coli - derived recombinant Human POLR3G (Met1-Met110).
Accession	O15318
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-POLR3G polyclonal antibody at 1ug/mL on various samples.
Lane 1: K562 cell lysate