


Product Information	
Catalog No.	PA14649
Product name	Anti-Human FAM20B Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Glycosaminoglycan xylosylkinase, 2.7.1.-, Xylose kinase, FAM20B, KIAA0475
Immunogen	E. coli - derived recombinant Human FAM20B (Gln61-Leu409).
Accession	O75063
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-FAM20B polyclonal antibody at FAM20B on various samples.

Lane 1: Mouse brain lysate

Lane 2: Mouse cerebellum lysate

Lane 3: Rat brain lysate

Lane 4: Rat cerebellum lysate