

Product Information	
Catalog No.	PA14656
Product name	Anti-Human FKBP6 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Inactive peptidyl-prolyl cis-trans isomerase FKBP6, Inactive PPIase FKBP6, 36 kDa FK506-binding protein, 36 kDa FKBP, FKBP-36, FK506-binding protein 6, FKBP-6, Immunophilin FKBP36, FKBP6, FKBP36
Immunogen	E. coli - derived recombinant Human FKBP6 (Ser20-Ser327).
Accession	O75344
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.