

Product Information	
Catalog No.	PA14666
Product name	Anti-Human TAF5L Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L, TAF5L, PCAF-associated factor 65 beta, PAF65-beta, TAF5L, PAF65B
Immunogen	E. coli - derived recombinant Human TAF5L (Thr257-His589).
Accession	O75529
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.