


Product Information	
Catalog No.	PA14675
Product name	Anti-Human SLC25A12 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Electrogenic aspartate/glutamate antiporter SLC25A12, mitochondrial, Araceli hiperlarga, Aralar, Aralar1, Mitochondrial aspartate glutamate carrier 1, Solute carrier family 25 member 12, SLC25A12, AGC1, ARALAR1
Immunogen	E. coli - derived recombinant Human SLC25A12 (Gly11-Gln311).
Accession	O75746
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-SLC25A12 polyclonal antibody at SLC25A12 on various samples.

- Lane 1: Jurkat cell lysate
- Lane 2: A431 cell lysate
- Lane 3: Mouse brain lysate
- Lane 4: Mouse skeletal muscle lysate
- Lane 5: Rat brain lysate
- Lane 6: Rat skeletal muscle lysate