


Product Information	
Catalog No.	PA14707
Product name	Anti-Human FKBP9 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Peptidyl-prolyl cis-trans isomerase FKBP9, PPIase FKBP9, 5.2.1.8, 63 kDa FK506-binding protein, 63 kDa FKBP, FKBP-63, FK506-binding protein 9, FKBP-9, Rotamase, FKBP9, FKBP60, FKBP63
Immunogen	E. coli - derived recombinant Human FKBP9 (Ala25-Leu570).
Accession	O95302
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-FKBP9 polyclonal antibody at FKBP9 on various samples.
Lane 1: NIH3T3 cell lysate
Lane 2: LO2 cell lysate
Lane 3: Mouse liver lysate