


Product Information	
Catalog No.	PA14718
Product name	Anti-Human ETHE1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Persulfide dioxygenase ETHE1, mitochondrial, 1.13.11.18, Ethylmalonic encephalopathy protein 1, Hepatoma subtracted clone one protein, Sulfur dioxygenase ETHE1, ETHE1, HSCO
Immunogen	E. coli - derived recombinant Human ETHE1 (Ile23-Ala254).
Accession	O95571
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-ETHE1 polyclonal antibody at ETHE1 on various samples.

Lane 1: Mouse kidney lysate

Lane 2: Mouse colon lysate

Lane 3: Rat liver lysate

Lane 4: Rat kidney lysate