


Product Information	
Catalog No.	PA14723
Product name	Anti-Human ESRRB Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	ERR-beta, NR3B2, Estrogen-related receptor beta, ESRL2, Estrogen receptor-like 2, ERR beta-2, ESRRB, Nuclear receptor subfamily 3 group B member 2, Steroid hormone receptor ERR2, ERRB2
Immunogen	E. coli - derived recombinant Human ESRRB (Gly235-Val433).
Accession	O95718
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-ESRRB polyclonal antibody at 1ug/mL on various samples.
Lane 1: Mouse heart lysate
Lane 2: Rat heart lysate