

Product Information	
Catalog No.	PA14747
Product name	Anti-Human PRH1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Salivary acidic proline-rich phosphoprotein 1/2, Db-s, PRP-1/PRP-2, Parotid acidic protein, Pa, Parotid double-band protein, Parotid isoelectric focusing variant protein, PIF-S, Parotid proline-rich protein 1/2, Pr1/Pr2, Protein C, Salivary acidic proline-rich phosphoprotein 1/2, Salivary acidic proline-rich phosphoprotein 3/4, Db-F, PIF-F, PRP-3/PRP-4, Protein A, Peptide P-C, PRH1, PRH2
Immunogen	E. coli - derived recombinant Human PRH1 (Gln17-Gln166).
Accession	P02810
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.