


Product Information	
Catalog No.	PA14784
Product name	Anti-Human HAPLN1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Hyaluronan and proteoglycan link protein 1, Cartilage-linking protein 1, Cartilage-link protein, Proteoglycan link protein, HAPLN1, CRTL1
Immunogen	E. coli - derived recombinant Human HAPLN1 (Gly37-Asn354).
Accession	P10915
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-HAPLN1 polyclonal antibody at 1ug/mL on various samples.

Lane 1: HepG2 cell lysate
Lane 2: Caco2 cell lysate
Lane 3: HUVEC cell lysate
Lane 4: SW480 cell lysate
Lane 5: U2OS cell lysate