


Product Information	
Catalog No.	PA14796
Product name	Anti-Human GTF2F2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	General transcription factor IIF subunit 2, General transcription factor IIF 30 kDa subunit, Transcription initiation factor IIF subunit beta, TFIIIF-beta, Transcription initiation factor RAP30, GTF2F2, RAP30
Immunogen	E. coli - derived recombinant Human GTF2F2 (Met1-Asp249).
Accession	P13984
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-GTF2F2 polyclonal antibody at 1ug/mL on various samples.

Lane 1: Hela cell lysate
Lane 2: K562 cell lysate
Lane 3: PC-3 cell lysate
Lane 4: 293T cell lysate