

Product Information	
Catalog No.	PA14853
Product name	Anti-Human EXOSC9 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Exosome complex component RRP45, Autoantigen PM/Scl 1, Exosome component 9, P75 polymyositis-scleroderma overlap syndrome-associated autoantigen, Polymyositis/scleroderma autoantigen 1, Polymyositis/scleroderma autoantigen 75 kDa, PM/Scl-75, EXOSC9, PMSCL1
Immunogen	E. coli - derived recombinant Human EXOSC9 (Met106-Ala295).
Accession	Q06265
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.