


Product Information	
Catalog No.	PA14942
Product name	Anti-Human HEPACAM2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	HEPACAM family member 2, Mitotic kinetics regulator, HEPACAM2, MIK1, UNQ305/PRO346
Immunogen	E. coli - derived recombinant Human HEPACAM2 (Leu32-Pro351).
Accession	A8MVW5
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-HEPACAM2 polyclonal antibody at 1ug/mL on various samples.
Lane 1: NCI-H460 cell lysate

Western blot