


Product Information	
Catalog No.	PA14957
Product name	Anti-Human SUGCT Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Succinate--hydroxymethylglutarate CoA-transferase, 2.8.3.13, Dermal papilla-derived protein 13, SuccinylCoA:glutarate-CoA transferase, SUGCT, C7orf10, DERP13
Immunogen	E. coli - derived recombinant Human SUGCT (Gln32-His438).
Accession	Q9HAC7
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-SYGCT polyclonal antibody at 1ug/mL on various samples.
Lane 1: Mouse kidney lysate

Western blot