


Product Information	
Catalog No.	PA14975
Product name	Anti-Mouse CD172a/SIRPA Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Tyrosine-protein phosphatase non-receptor type substrate 1, SHP substrate 1, SHPS-1, Brain Ig-like molecule with tyrosine-based activation motifs, Bit, CD172 antigen-like family member A, Inhibitory receptor SHPS-1, MyD-1 antigen, Signal-regulatory protein alpha-1, Sirp-alpha-1, mSIRP-alpha1, p84, CD172a, Sirpa, Bit, Myd1, Ptpns1, Shps1, Sirp
Immunogen	E. coli - derived recombinant Mouse CD172a/SIRPA (Gly31-Gly350).
Accession	P97797
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-SIRPA polyclonal antibody at SIRPA on various samples.
Lane 1: Mouse brain lysate
Lane 2: Mouse cerebellum lysate
Lane 3: Rat brain lysate