


Product Information	
Catalog No.	PA14994
Product name	Anti-Mouse ADA/Adenosine deaminase Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Adenosine deaminase, 3.5.4.4, Adenosine aminohydrolase, Ada
Immunogen	E. coli - derived recombinant Mouse ADA/Adenosine deaminase (Met1-Gln352).
Accession	P03958
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-ADA polyclonal antibody at ADA on various samples.

- Lane 1: Jurkat cell lysate
- Lane 2: MOLT-4 cell lysate
- Lane 3: Mouse thymus lysate
- Lane 4: Mouse stomach lysate
- Lane 5: Rat spleen lysate