

Product Information	
Catalog No.	PA15107
Product name	Anti-ASFV CP2475L (p37) Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	African swine fever virus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Polyprotein pp220, 220 kDa polyprotein, p34, p14, p37, p150, p5, Ba71V-92, CP2475L
Immunogen	E. coli - derived recombinant ASFV CP2475L (p37) (Ala523-Gly893).
Accession	Q08358
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.