

Product Information	
Catalog No.	PA15142
Product name	Anti-BCCV Gn/Glycoprotein N Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Black Creek Canal orthohantavirus (BCCV) (Black Creek Canal virus)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Envelopment polyprotein, Glycoprotein precursor, M polyprotein, Glycoprotein N, Gn, Glycoprotein G1, Glycoprotein C, Gc, Glycoprotein G2, GP
Immunogen	E. coli - derived recombinant BCCV Gn/Glycoprotein N (Phe19-Val482).
Accession	P0DTJ0
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.