

Product Information	
Catalog No.	PA15150
Product name	Anti-HAdV-5 E1B-55K Polyclonal Antibody
Applications	ELISA, WB
Species reactivity	Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	E1B 55 kDa protein, E1B-55K, E1B protein, large T-antigen, E1B-495R
Immunogen	E. coli - derived recombinant HAdV-5 E1B-55K (Glu2-Gln79).
Accession	P03243
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.