

Product Information	
Catalog No.	PA15155
Product name	Anti-HAdV-4 Protease/L3 Polyclonal Antibody
Applications	ELISA, WB
Species reactivity	Human adenovirus E serotype 4 (HAdV-4) (Human adenovirus 4)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Protease, 3.4.22.39, Adenain, Adenovirus protease, AVP, Adenovirus proteinase, Endoprotease, L3
Immunogen	E. coli - derived recombinant HAdV-4 Protease/L3 (Met1-Met201).
Accession	P07885
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.